


HA Tag Antibody [clone BC-8] (RQ4678)

Catalog No.	Formulation	Size
RQ4678	Antibody in PBS with 0.02% sodium azide, 50% glycerol and 0.4-0.5mg/ml BSA	100 ul

[Bulk quote request](#)

Availability	1-2 weeks
Species Reactivity	HA Tag
Format	Purified
Host	Rabbit
Clonality	Rabbit Monoclonal
Isotype	Rabbit IgG
Clone Name	BC-8
Purity	Affinity purified
Applications	Western Blot : 1:3000-1:10000
Limitations	This HA Tag antibody is available for research use only.



Western blot testing of an HA-tagged protein with HA Tag antibody.

Description

Epitope tags consisting of short sequences recognized by well-characterized monoclonal antibodies have been widely used in the study of protein expression in various systems. NSJ's anti-HA monoclonal antibody (Clone BC-8) provides a simple solution to detect the expression of an HA-tagged protein in cells.

Application Notes

Optimal dilution of the HA Tag antibody should be determined by the researcher.

Immunogen

A synthetic peptide specific to Human influenza hemagglutinin (HA) was used as the immunogen for the HA Tag antibody.

Storage

Store the HA Tag antibody at -20oC.