


## HAO2 Antibody / Hydroxyacid oxidase 2 (R33671)

Catalog No.	Formulation	Size
R33671-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

<b>Availability</b>	1-3 business days
<b>Species Reactivity</b>	Human
<b>Format</b>	Antigen affinity purified
<b>Host</b>	Goat
<b>Clonality</b>	Polyclonal (goat origin)
<b>Isotype</b>	Goat Ig
<b>Purity</b>	Antigen affinity
<b>Gene ID</b>	51179
<b>Applications</b>	Western Blot : 0.3-1ug/ml Enzyme ImmunoAssay : binds recombinant protein (use at 5-10ug/ml) ELISA (peptide) LOD : 1:64000
<b>Limitations</b>	This HAO2 antibody is available for research use only.



Western blot testing of human kidney lysate with HAO2 antibody at 0.3ug/ml. Predicted molecular weight: ~39 kDa.

### Description

Additional name(s) for this target protein: 2 hydroxy acid oxidase; HAOX2

### Application Notes

Optimal dilution of the HAO2 antibody should be determined by the researcher.

## **Immunogen**

Amino acids DVSEVDTRTTIQGE were used as the immunogen for this HAO2 antibody.

## **Storage**

Aliquot and store the HAO2 antibody at -20oC.