


HAO2 Antibody / Hydroxyacid oxidase 2 (R33671)

Catalog No.	Formulation	Size
R33671-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide $$	100 ug

Bulk quote request

Availability	1-3 business days
Species Reactivity	Human
Format	Antigen affinity purified
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	51179
Applications	Western Blot : 0.3-1ug/ml Enzyme ImmunoAssay : binds recombinant protein (use at 5-10ug/ml) ELISA (peptide) LOD : 1:64000
Limitations	This HAO2 antibody is available for research use only.



Description

Additional name(s) for this target protein: 2 hydroxy acid oxidase; HAOX2

Application Notes

Optimal dilution of the HAO2 antibody should be determined by the researcher.

Immunogen

Amino acids DVSEVDTRTTIQGE were used as the immunogen for this HAO2 antibody. **Storage** Aliquot and store the HAO2 antibody at -20oC.