


HAVCR2 Antibody / TIM-3 [clone 1827CT151.74.68] (F54121)

Catalog No.	Formulation	Size
F54121-0.2ML	In 1X PBS, pH 7.4, with 0.09% sodium azide	0.2 ml
F54121-0.05ML	In 1X PBS, pH 7.4, with 0.09% sodium azide	0.05 ml

Bulk quote request

Availability	1-3 business days
Species Reactivity	Human
Format	Purified
Host	Mouse
Clonality	Monoclonal (mouse origin)
Isotype	Mouse IgG2b, kappa
Clone Name	1827CT151.74.68
Purity	Protein G purified
UniProt	Q8TDQ0
Applications	Western Blot : 1:4000
Limitations	This HAVCR2 antibody is available for research use only.



Western blot testing of HAVCR2 antibody at 1:4000: Lane 1) human HepG2, 2) Daudi and 3) kidney lysate. Expected molecular weight: 33-70 kDa depending on glycosylation level.

Description

Regulates macrophage activation. Inhibits T-helper type 1 lymphocyte (Th1)-mediated auto- and alloimmune responses and promotes immunological tolerance. May be also involved in T-cell homing. Receptor for LGALS9. [UniProt]

Application Notes

The stated application concentrations are suggested starting points. Titration of the HAVCR2 antibody may be required

due to differences in protocols and secondary/substrate sensitivity.

Immunogen

Recombinant human HAVCR2 was used as the immunogen for the HAVCR2 antibody.

Storage

Aliquot the HAVCR2 antibody and store frozen at -20oC or colder. Avoid repeated freeze-thaw cycles.