


## HCN3 Antibody (R36027)

Catalog No.	Formulation	Size
R36027-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human, Mouse
Predicted Reactivity	Cow, Dog, Pig, Rat
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	57657
Applications	Western Blot : 0.01-0.03ug/ml ELISA (peptide) LOD : 1:8000
Limitations	This HCN3 antibody is available for research use only.



Western blot testing of mouse brain lysate with HCN3 antibody at 0.01ug/ml. Predicted molecular weight: ~86kDa.

## Description

Additional name(s) for this target protein: Hyperpolarization activated cyclic nucleotide-gated potassium channel 3

## Application Notes

Optimal dilution of the HCN3 antibody should be determined by the researcher.

## **Immunogen**

Amino acids SQPSLPQRATGDGS were used as the immunogen for this HCN3 antibody.

## **Storage**

Aliquot and store the HCN3 antibody at -20oC.