


## HDAC3 Antibody (R34813)

Catalog No.	Formulation	Size
R34813-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Predicted Reactivity	Cow, Dog, Mouse, Pig, Rat
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	8841
Applications	Western Blot : 0.3-1ug/ml ELISA (peptide) LOD : 1:128000
Limitations	This HDAC3 antibody is available for research use only.



Western blot testing of Jurkat nuclear lysate with HDAC3 antibody at 0.1ug/ml. Predicted molecular weight ~49kDa.

### Description

Additional name(s) for this target protein: Histone deacetylase 3

### Application Notes

Optimal dilution of the HDAC3 antibody should be determined by the researcher.

## **Immunogen**

Amino acids NHAPSVQIHDVPAD were used as the immunogen for this HDAC3antibody.

## **Storage**

Aliquot and store the HDAC3 antibody at -20oC.