


HDAC6 Antibody (R33139)

Catalog No.	Formulation	Size
R33139-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human, Mouse
Predicted Reactivity	Cow, Dog, Rat
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	10013
Applications	Western Blot : 1-3ug/ml ELISA (peptide) LOD : 1:32000
Limitations	This HDAC6 antibody is available for research use only.



Western blot testing of HeLa nuclear lysate with HDAC6 antibody at 1ug/ml. Expected molecular weight 130~160kDa.

Description

Additional name(s) for this target protein: Histone deacetylase 6

Application Notes

Optimal dilution of the HDAC6 antibody should be determined by the researcher.

Immunogen

Amino acids AHQNKFGEDMPHPH were used as the immunogen for this HDAC6 antibody.

Storage

Aliquot and store the HDAC6 antibody at -20oC.