


HDGFRP3 Antibody (R33898)

Catalog No.	Formulation	Size
R33898-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Predicted Reactivity	Cow, Dog
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	50810
Applications	Western Blot : 0.3-1ug/ml ELISA (peptide) LOD : 1:128000
Limitations	This HDGFRP3 antibody is available for research use only.



Western blot testing of human amygdala lysate with HDGFRP3 antibody at 0.3ug/ml.
Predicted molecular weight: ~23kDa.

Description

Additional name(s) for this target protein: Hepatoma-derived growth factor-related protein 3

Application Notes

Optimal dilution of the HDGFRP3 antibody should be determined by the researcher.

Immunogen

Amino acids ENNPGVKFTGYQAIQ were used as the immunogen for this HDGFRP3 antibody.

Storage

Aliquot and store the HDGFRP3 antibody at -20oC.