


# **HDGFRP3 Antibody (R33898)**

Catalog No.	Formulation	Size
R33898-100UG	0.5~mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

#### **Bulk quote request**

Availability	1-3 business days
Species Reactivity	Human
Predicted Reactivity	Cow, Dog
Format	Antigen affinity purified
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	50810
Applications	Western Blot : 0.3-1ug/ml ELISA (peptide) LOD : 1:128000
Limitations	This HDGFRP3 antibody is available for research use only.



## **Description**

Additional name(s) for this target protein: Hepatoma-derived growth factor-related protein 3

### **Application Notes**

Optimal dilution of the HDGFRP3 antibody should be determined by the researcher.

### **Immunogen**

Amino acids ENNPGVKFTGYQAIQ were used as the immunogen for this HDGFRP3 antibody. **Storage** Aliquot and store the HDGFRP3 antibody at -20oC.