


HES4 Antibody (R35757)

Catalog No.	Formulation	Size
R35757-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	57801
Applications	Western Blot : 1-3ug/ml ELISA (peptide) LOD : 1:128000
Limitations	This HES4 antibody is available for research use only.



Western blot testing of HepG2 nuclear lysate with HES4 antibody at 1ug/ml. Predicted molecular weight: ~26kDa.

Description

Additional name(s) for this target protein: Hairy and enhancer of split 4

Application Notes

Optimal dilution of the HES4 antibody should be determined by the researcher.

1. This HES4 antibody is specific for isoform 1.

Immunogen

Amino acids RSAAEHRKVGSRP were used as the immunogen for this HES4 antibody.

Storage

Aliquot and store the HES4 antibody at -20oC.