


HES7 Antibody (F51870)

Catalog No.	Formulation	Size
F51870-0.4ML	In 1X PBS, pH 7.4, with 0.09% sodium azide	0.4 ml
F51870-0.08ML	In 1X PBS, pH 7.4, with 0.09% sodium azide	0.08 ml

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human, Mouse
Format	Antigen affinity purified
Host	Rabbit
Clonality	Polyclonal (rabbit origin)
Isotype	Rabbit Ig
Purity	Antigen affinity
UniProt	Q9BYE0
Applications	Western Blot : 1:1000
Limitations	This HES7 antibody is available for research use only.



Description

HES7 is a member of the hairy and enhancer of split family of bHLH transcription factors. The protein functions as a transcriptional repressor, and is implicated in correct patterning of the axial skeleton.

Application Notes

Titration of the HES7 antibody may be required due to differences in protocols and secondary/substrate sensitivity.

Immunogen

A portion of amino acids 26-55 from the human protein was used as the immunogen for this HES7 antibody.

Storage

Aliquot the HES7 antibody and store frozen at -20oC or colder. Avoid repeated freeze-thaw cycles.