


Histone H1.10 Antibody / H1-10 (RQ8788)


Catalog No.	Formulation	Size
RQ8788	0.5mg/ml if reconstituted with 0.2ml sterile DI water	100 ug

[Bulk quote request](#)

Availability	1-3 days
Species Reactivity	Human
Format	Antigen affinity purified
Host	Rabbit
Clonality	Polyclonal (rabbit origin)
Isotype	Rabbit IgG
Purity	Antigen affinity chromatography
Buffer	Lyophilized from 1X PBS with 2% Trehalose
UniProt	Q92522
Localization	Nuclear
Applications	Western Blot : 1-2ug/ml Immunofluorescence : 5ug/ml ELISA : 0.1-0.5ug/ml
Limitations	This Histone H1.10 antibody is available for research use only.



Immunofluorescent staining of FFPE human HeLa cells with Histone H1.10 antibody (red) and Beta Tubulin mAb (green). HIER: steam section in pH6 citrate buffer for 20 min.



Western blot testing of human 1) HeLa, 2) 293T and 3) U-2 OS cell lysate with Histone H1.10 antibody. Predicted molecular weight ~22 kDa.

Description

Histone H1x is a protein that in humans is encoded by the H1FX gene. Histones are basic nuclear proteins that are responsible for the nucleosome structure of the chromosomal fiber in eukaryotes. Nucleosomes consist of approximately 146 bp of DNA wrapped around a histone octamer composed of pairs of each of the four core histones (H2A, H2B, H3, and H4). The chromatin fiber is further compacted through the interaction of a linker histone, H1, with the DNA between the nucleosomes to form higher order chromatin structures. This gene encodes a replication-independent histone that is a member of the histone H1 family.

Application Notes

Optimal dilution of the Histone H1.10 antibody should be determined by the researcher.

Immunogen

An E.coli-derived human recombinant protein (amino acids M1-K106) was used as the immunogen for the Histone H1.10 antibody.

Storage

After reconstitution, the Histone H1.10 antibody can be stored for up to one month at 4oC. For long-term, aliquot and store at -20oC. Avoid repeated freezing and thawing.