


HLA-DQA2 Antibody (R35277)

Catalog No.	Formulation	Size
R35277-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	3118; 3117
Applications	Western Blot : 1-3ug/ml ELISA (peptide) LOD : 1:32000
Limitations	This HLA-DQA2 antibody is available for research use only.



Western blot testing of human bone marrow lysate with HLA-DQA2 antibody at 1ug/ml.
Predicted molecular weight: ~28kDa.

Description

Additional name(s) for this target protein: Major histocompatibility complex, class II, DQ alpha 2

Application Notes

Optimal dilution of the HLA-DQA2 antibody should be determined by the researcher.

1. This HLA-DQA2 antibody will also detect the highly homologous HLA-DQA1.

Immunogen

Amino acids QGLRSVGASRH were used as the immunogen for this HLA-DQA2 antibody.(1)

Storage

Aliquot and store the HLA-DQA2 antibody at -20oC.