


HLTF Antibody (R34582)

Catalog No.	Formulation	Size
R34582-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Predicted Reactivity	Cow, Dog, Mouse, Rat
Format	Antigen affinity purified
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	6596
Applications	Western Blot : 0.5-1.5ug/ml ELISA (peptide) LOD : 1:16000
Limitations	This HLTF antibody is available for research use only.



Western blot testing of Jurkat lysate with HLTF antibody at 0.5ug/ml. Predicted molecular weight: ~114kDa.

Description

Additional name(s) for this target protein: Helicase-like transcription factor, SMARCA3

Application Notes

Optimal dilution of the HLTF antibody should be determined by the researcher.

Immunogen

Amino acids KQAKINEIRTLIDL were used as the immunogen for this HLTF antibody.

Storage

Aliquot and store the HLTF antibody at -20oC.