


HMBS Antibody (R34224)

Catalog No.	Formulation	Size
R34224-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	3145
Applications	Western Blot : 0.1-0.3ug/ml ELISA (peptide) LOD : 1:32000
Limitations	This HMBS antibody is available for research use only.



Western blot testing of K562 lysate with HMBS antibody at 0.1ug/ml. Predicted molecular weight: ~39 kDa.

Description

Additional name(s) for this target protein: Hydroxymethylbilane synthase

Application Notes

Optimal dilution of the HMBS antibody should be determined by the researcher.

Immunogen

Amino acids HLEFRSIRGNLNTR were used as the immunogen for this HMBS antibody.

Storage

Aliquot and store the HMBS antibody at -20oC.