


HMG4 Antibody / HMGB3 (R34430)

Catalog No.	Formulation	Size
R34430-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Predicted Reactivity	Cow, Mouse
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	3149
Applications	Western Blot : 0.01-0.03ug/ml ELISA (peptide) LOD : 1:64000
Limitations	This HMG4 antibody is available for research use only.



Western blot testing of human lung lysate with HMG4 antibody at 0.01ug/ml. Predicted molecular weight: ~23kDa.

Description

Additional name(s) for this target protein: High-mobility group (nonhistone chromosomal) protein 4, High-mobility group box 3

Application Notes

Optimal dilution of the HMG4 antibody should be determined by the researcher.

Immunogen

Amino acids KFDGAKGPAKVARKK were used as the immunogen for this HMG4 antibody.

Storage

Aliquot and store the HMG4 antibody at -20oC.