


## HMGA2 Antibody (R34342)

Catalog No.	Formulation	Size
R34342-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

<b>Availability</b>	1-3 business days
<b>Species Reactivity</b>	Human
<b>Predicted Reactivity</b>	Dog
<b>Format</b>	Antigen affinity purified
<b>Host</b>	Goat
<b>Clonality</b>	Polyclonal (goat origin)
<b>Isotype</b>	Goat Ig
<b>Purity</b>	Antigen affinity
<b>Gene ID</b>	8091
<b>Applications</b>	Western Blot : 1-3ug/ml ELISA (peptide) LOD : 1:8000
<b>Limitations</b>	This HMGA2 antibody is available for research use only.



Western blot testing of human heart lysate with HMGA2 antibody at 1ug/ml. Predicted molecular weight ~12kDa but routinely observed at ~18kDa.

### Description

Additional name(s) for this target protein: High mobility group AT-hook 2, BABL, HMGI-C

### Application Notes

Optimal dilution of the HMGA2 antibody should be determined by the researcher.

## **Immunogen**

Amino acids KAAQKKAETGEK were used as the immunogen for this HMGA2 antibody.

## **Storage**

Aliquot and store the HMGA2 antibody at -20oC.