


HMGCs1 Antibody (R36420)

Catalog No.	Formulation	Size
R36420-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	3157
Applications	Western Blot : 1-3ug/ml ELISA (peptide) LOD : 1:128000
Limitations	This HMGCs1 antibody is available for research use only.



Western blot testing of A) HeLa, B) HepG2 and C) K562 cell lysate using HMGCs1 antibody at 1ug/ml. Predicted/observed molecular weight: ~57kDa.

Description

Additional name(s) for this target protein: 3-hydroxy-3-methylglutaryl-CoA synthase 1

Application Notes

Optimal dilution of the HMGCs1 antibody should be determined by the researcher.

Immunogen

Amino acids KASSELSQKTK were used as the immunogen for this HMGCS1 antibody.

Storage

Aliquot and store the HMGCS1 antibody at -20oC.