


hnRNP D Antibody / HNRNPD / AUF1 [clone 4F3] (RQ5859)


Catalog No.	Formulation	Size
RQ5859	0.5mg/ml if reconstituted with 0.2ml sterile DI water	100 ug

Bulk quote request


Availability	1-3 business days
Species Reactivity	Human, Mouse, Rat
Format	Antigen affinity purified
Host	Mouse
Clonality	Monoclonal (mouse origin)
Isotype	Mouse IgG1
Clone Name	4F3
Purity	Affinity purified
Buffer	Lyophilized from 1X PBS with 2% Trehalose and 0.025% sodium azide
UniProt	Q14103
Localization	Nuclear, cytoplasmic
Applications	Western Blot : 0.5-1ug/ml Immunohistochemistry : 1-2ug/ml Immunofluorescence : 2-4ug/ml Immunocytochemistry : 2-4ug/ml Flow Cytometry : 1-3ug/million cells
Limitations	This hnRNP D antibody is available for research use only.




IHC staining of FFPE human rectal cancer with hnRNP D antibody. HIER: boil tissue sections in pH8 EDTA for 20 min and allow to cool before testing.




Immunofluorescent staining of FFPE human rectal cancer with hnRNP D antibody. HIER: boil tissue sections in pH8 EDTA for 20 min and allow to cool before testing.




IHC staining of FFPE rat brain with hnRNP D antibody. HIER: boil tissue sections in pH8 EDTA for 20 min and allow to cool before testing.



Western blot testing of human 1) HEK293, 2) K562, 3) A431, 4) HepG2, 5) placenta, 6) rat brain and 7) mouse brain lysate with hnRNP D antibody. Expected molecular weight: multiple bands from 37-45 kDa.



Flow cytometry testing of human SiHa cells with hnRNP D antibody at 1ug/million cells (blocked with goat sera); Red=cells alone, Green=isotype control, Blue= hnRNP D antibody.



Flow cytometry testing of human U937 cells with hnRNP D antibody at 1ug/million cells (blocked with goat sera); Red=cells alone, Green=isotype control, Blue= hnRNP D antibody.

Description

Heterogeneous nuclear ribonucleoprotein D0 (HNRNPD) also known as AU-rich element RNA-binding protein 1 (AUF1) is a protein that in humans is encoded by the HNRNPD gene. It is mapped to 4q21.22. This gene belongs to the subfamily of ubiquitously expressed heterogeneous nuclear ribonucleoproteins (hnRNPs). The hnRNPs are nucleic acid binding proteins and they complex with heterogeneous nuclear RNA (hnRNA). These proteins are associated with pre-mRNAs in the nucleus and appear to influence pre-mRNA processing and other aspects of mRNA metabolism and transport. While all of the hnRNPs are present in the nucleus, some seem to shuttle between the nucleus and the cytoplasm. The hnRNP proteins have distinct nucleic acid binding properties. The protein encoded by this gene has two repeats of quasi-RRM domains that bind to RNAs. It localizes to both the nucleus and the cytoplasm. This protein is implicated in the regulation

of mRNA stability. Alternative splicing of this gene results in four transcript variants.

Application Notes

Optimal dilution of the hnRNP D antibody should be determined by the researcher.

Immunogen

Recombinant human protein (amino acids E88-N246) was used as the immunogen for the hnRNP D antibody.

Storage

After reconstitution, the hnRNP D antibody can be stored for up to one month at 4°C. For long-term, aliquot and store at -20°C. Avoid repeated freezing and thawing.