


Homeobox protein cut-like 1 Antibody / CUX1 (RQ7069)


Catalog No.	Formulation	Size
RQ7069	0.5mg/ml if reconstituted with 0.2ml sterile DI water	100 ug

Bulk quote request


Availability	1-3 business days
Species Reactivity	Human, Mouse
Format	Antigen affinity purified
Host	Rabbit
Clonality	Polyclonal (rabbit origin)
Isotype	Rabbit IgG
Purity	Antigen affinity purified
Buffer	Lyophilized from 1X PBS with 2% Trehalose
UniProt	P39880
Localization	Nuclear, cytoplasmic (Golgi)
Applications	Western Blot : 0.5-1ug/ml Immunohistochemistry (FFPE) : 2-5ug/ml Flow Cytometry : 1-3ug/million cells Direct ELISA : 0.1-0.5ug/ml
Limitations	This Homeobox protein cut-like 1 antibody is available for research use only.




IHC staining of FFPE human laryngeal squamous cell carcinoma tissue with Homeobox protein cut-like 1 antibody. HIER: boil tissue sections in pH8 EDTA for 20 min and allow to cool before testing.



IHC staining of FFPE human thyroiditis tissue with Homeobox protein cut-like 1 antibody.
HIER: boil tissue sections in pH8 EDTA for 20 min and allow to cool before testing.




IHC staining of FFPE human placental tissue with Homeobox protein cut-like 1 antibody.
HIER: boil tissue sections in pH8 EDTA for 20 min and allow to cool before testing.



IHC staining of FFPE mouse colon tissue with Homeobox protein cut-like 1 antibody.
HIER: boil tissue sections in pH8 EDTA for 20 min and allow to cool before testing.



Western blot testing of 1) human K562, 2) human HeLa, 3) human MCF7, 4) human SH-SY5Y and 5) mouse stomach tissue lysate with Homeobox protein cut-like 1 antibody.
Expected molecular weight ~200 kDa with smaller isoforms.



Flow cytometry testing of human HepG2 cells with Homeobox protein cut-like 1 antibody at 1 μg/million cells (blocked with goat sera); Red=cells alone, Green=isotype control, Blue= Homeobox protein cut-like 1 antibody.

Description

The protein encoded by this gene is a member of the homeodomain family of DNA binding proteins. It may regulate gene expression, morphogenesis, and differentiation and it may also play a role in the cell cycle progression. Several alternatively spliced transcript variants encoding different isoforms have been identified.

Application Notes

Optimal dilution of the Homeobox protein cut-like 1 antibody should be determined by the researcher.

Immunogen

Recombinant human protein (amino acids M1-E220) was used as the immunogen for the CUX1 antibody.

Storage

After reconstitution, the Homeobox protein cut-like 1 antibody can be stored for up to one month at 4oC. For long-term, aliquot and store at -20oC. Avoid repeated freezing and thawing.