


HOXA1 Antibody (Center Region) (R35385)

Catalog No.	Formulation	Size
R35385-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Predicted Reactivity	Mouse, Rat, Dog, Cow
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	3198
Applications	Western Blot : 1-3ug/ml ELISA (peptide) LOD : 1:8000
Limitations	This HOXA1 antibody is available for research use only.



HEK293 lysate overexpressing human HOXA1-MYC probed with HOXA1 [1ug/ml] in Lane A and anti-MYC tag [1/1000] in lane C. Mock-transfected HEK293 lysate probed with HOXA1 antibody [1ug/ml] in Lane B. Predicted molecular weight: ~37/15kDa (isoforms a/b).

Description

Additional name(s) for this target protein: Homeobox A1

Application Notes

Optimal dilution of the HOXA1 antibody should be determined by the researcher.

1. This HOXA1 antibody is specific for isoforms a and b.

Immunogen

Amino acids QPATYQTSGN were used as the immunogen for this HOXA1 antibody.

Storage

Aliquot and store the HOXA1 antibody at -20oC.