


HOXA4 Antibody (R30682)


Catalog No.	Formulation	Size
R30682	0.5mg/ml if reconstituted with 0.2ml sterile DI water	100 ug

[Bulk quote request](#)


Availability	1-3 business days
Species Reactivity	Human, Mouse, Rat
Format	Antigen affinity purified
Host	Rabbit
Clonality	Polyclonal (rabbit origin)
Isotype	Rabbit IgG
Purity	Antigen affinity
Buffer	Lyophilized from 1X PBS with 2.5% BSA and 0.025% sodium azide/thimerosal
UniProt	Q00056
Applications	Western Blot : 0.5-1ug/ml IHC (FFPE) : 0.5-1ug/ml Immunocytochemistry : 0.5-1ug/ml
Limitations	This HOXA4 antibody is available for research use only.



Western blot testing of HOXA4 antibody and Lane 1: SW620; 2: SW620; 3: PC-12 cell nuclear lysate. Predicted molecular weight: ~35kDa.



IHC-P: HOXA4 antibody testing of human breast cancer tissue. HIER: steamed with pH6 citrate buffer.



ICC testing of HeLa cells with HOXA4 antibody.

Description

HOXA4 (HOMEODOMAIN A4), also known as HOX1D or HOMOLOG OF, is a protein that in humans is encoded by the HOXA4 gene, which is also part of the A cluster on chromosome 7. HOXA4 encodes a DNA-binding transcription factor which may regulate gene expression, morphogenesis, and differentiation. Its cytogenetic location is 7q15.2. In vertebrates, the genes encoding the class of transcription factors called homeobox genes are found in clusters named A, B, C, and D on four separate chromosomes. Expression of these proteins is spatially and temporally regulated during embryonic development.

Application Notes

The stated application concentrations are suggested starting amounts. Titration of the HOXA4 antibody may be required due to differences in protocols and secondary/substrate sensitivity.

Immunogen

An amino acid sequence from the N-terminus of human HOXA4 (NSNYIEPKFPPFEEYAQH) was used as the immunogen for this HOXA4 antibody.

Storage

After reconstitution, the HOXA4 antibody can be stored for up to one month at 4°C. For long-term, aliquot and store at -20°C. Avoid repeated freezing and thawing.