


HOXA4 Antibody (R35194)

Catalog No.	Formulation	Size
R35194-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Predicted Reactivity	Dog
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	3201
Applications	Western Blot : 0.1-0.3ug/ml ELISA (peptide) LOD : 1:32000
Limitations	This HOXA4 antibody is available for research use only.



Western blot testing of human colon lysate with HOXA4 antibody at 0.1ug/ml. Predicted molecular weight: ~35kDa.

Description

Additional name(s) for this target protein: Homeobox A4

Application Notes

Optimal dilution of the HOXA4 antibody should be determined by the researcher.

Immunogen

Amino acids QAKGPAHGLHASH were used as the immunogen for this HOXA4 antibody.

Storage

Aliquot and store the HOXA4 antibody at -20oC.