


# HOXA4 Antibody (R35194)

Catalog No.	Formulation	Size
R35194-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide $$	100 ug

#### **Bulk quote request**

Availability	1-3 business days
Species Reactivity	Human
Predicted Reactivity	Dog
Format	Antigen affinity purified
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	3201
Applications	Western Blot : 0.1-0.3ug/ml ELISA (peptide) LOD : 1:32000
Limitations	This HOXA4 antibody is available for research use only.



### **Description**

Additional name(s) for this target protein: Homeobox A4

## **Application Notes**

Optimal dilution of the HOXA4 antibody should be determined by the researcher.

#### **Immunogen**

Amino acids QAKGPAHGLHASH were used as the immunogen for this HOXA4 antibody. **Storage** Aliquot and store the HOXA4 antibody at -20oC.