


## HOXA5 Antibody (R33926)


| Catalog No.  | Formulation                                                                   | Size   |
|--------------|-------------------------------------------------------------------------------|--------|
| R33926-100UG | 0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide | 100 ug |

[Bulk quote request](#)

|                    |                                                         |
|--------------------|---------------------------------------------------------|
| Availability       | 1-3 business days                                       |
| Species Reactivity | Human, Mouse, Rat, Pig                                  |
| Format             | Antigen affinity purified                               |
| Host               | Goat                                                    |
| Clonality          | Polyclonal (goat origin)                                |
| Isotype            | Goat Ig                                                 |
| Purity             | Antigen affinity                                        |
| Gene ID            | 3202                                                    |
| Applications       | Western Blot : 1-3ug/ml<br>ELISA (peptide) LOD : 1:8000 |
| Limitations        | This HOXA5 antibody is available for research use only. |



Western blot testing of human [A], mouse [B] and rat [C] kidney lysates with HOXA5 antibody at 1ug/ml. Predicted molecular weight: ~29 kDa.



Western blot of pig kidney lysate with HOXA5 antibody at 1ug/ml. Predicted molecular weight: ~29kDa.

## Description

Additional name(s) for this target protein: Homeobox A5

## Application Notes

Optimal dilution of the HOXA5 antibody should be determined by the researcher.

## Immunogen

Amino acids EPRYSQPATS were used as the immunogen for this HOXA5 antibody.

## Storage

Aliquot and store the HOXA5 antibody at -20oC.