

HOXA7 Antibody (Center Region) (F54098)

Catalog No.	Formulation	Size
F54098-0.2ML	In 1X PBS, pH 7.4, with 0.09% sodium azide	0.2 ml
F54098-0.05ML	In 1X PBS, pH 7.4, with 0.09% sodium azide	0.05 ml

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Format	Antigen affinity purified
Host	Rabbit
Clonality	Polyclonal (rabbit origin)
Isotype	Rabbit Ig
Purity	Antigen affinity
UniProt	P31268
Applications	Western Blot : 1:2000
Limitations	This HOXA7 antibody is available for research use only.

Western blot testing of HOXA7 antibody at 1:2000 + human skeletal muscle lysate.
Predicted molecular weight ~25 kDa.

Description

Sequence-specific transcription factor which is part of a developmental regulatory system that provides cells with specific positional identities on the anterior-posterior axis. [UniProt]

Application Notes

The stated application concentrations are suggested starting points. Titration of the HOXA7 antibody may be required due to differences in protocols and secondary/substrate sensitivity.

Immunogen

A portion of amino acids 82-114 from human HOXA7 was used as the immunogen for the HOXA7 antibody.

Storage

Aliquot the HOXA7 antibody and store frozen at -20°C or colder. Avoid repeated freeze-thaw cycles.

