


HOXA9 Antibody (R34969)

Catalog No.	Formulation	Size
R34969-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human, Mouse
Predicted Reactivity	Cow, Dog, Rat
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	3205
Applications	Western Blot : 2-4ug/ml ELISA (peptide) LOD : 1:8000
Limitations	This HOXA9 antibody is available for research use only.



Western blot testing of mouse spleen lysate with HOXA9 antibody at 2ug/ml. Predicted molecular weight ~30kDa, both bands are blocked by the immunizing peptide.

Description

Additional name(s) for this target protein: Homeobox A9

Application Notes

Optimal dilution of the HOXA9 antibody should be determined by the researcher.

Immunogen

Amino acids PDFSPCSFQSKA were used as the immunogen for this HOXA9 antibody.

Storage

Aliquot and store the HOXA9 antibody at -20oC.