


# HOXA9 Antibody (R35395)

Catalog No.	Formulation	Size
R35395-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide $$	100 ug

### **Bulk quote request**

Availability	1-3 business days
Species Reactivity	Human, Mouse, Rat
Format	Antigen affinity purified
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	3205
Applications	Western Blot : 1-3ug/ml ELISA (peptide) LOD : 1:16000
Limitations	This HOXA9 antibody is available for research use only.



## **Description**

Additional name(s) for this target protein: Homeobox A9

### **Application Notes**

Optimal dilution of the HOXA9 antibody should be determined by the researcher.

### Immunogen

Amino acids QPPRQAATLAEHPD were used as the immunogen for this HOXA9 antibody.

# **Storage** Aliquot and store the HOXA9 antibody at -20oC.