

HOXB13 Antibody (R34502)

Catalog No.	Formulation	Size
R34502-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	10481
Applications	Western Blot : 0.1-0.3ug/ml ELISA (peptide) LOD : 1:128000
Limitations	This HOXB13 antibody is available for research use only.



Western blot testing of HeLa cell lysate with HOXB13 antibody at 0.1ug/ml. Predicted molecular weight: ~31kDa.

Description

Additional name(s) for this target protein: Homeobox B13

Application Notes

Optimal dilution of the HOXB13 antibody should be determined by the researcher.

Immunogen

Amino acids KKVLAQVKNSATP were used as the immunogen for this HOXB13 antibody.

Storage

Aliquot and store the HOXB13 antibody at -20oC.