


HOXB6 Antibody (R33923)


Catalog No.	Formulation	Size
R33923-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Predicted Reactivity	Cow, Dog, Mouse, Pig, Rat
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	3216
Applications	Western Blot : 1-3ug/ml Immunohistochemistry (FFPE) : 5-10ug/ml ELISA (peptide) LOD : 1:1000
Limitations	This HOXB6 antibody is available for research use only.



Western blot testing of HEK293 lysate with HOXB6 antibody at 1ug/ml. Predicted molecular weight: ~25kDa.



IHC staining of FFPE human lung with HOXB6 antibody at 5ug/ml. HIER: steamed with pH6 citrate buffer, AP-staining.

Description

Additional name(s) for this target protein: Homeobox B6

Application Notes

Optimal dilution of the HOXB6 antibody should be determined by the researcher.

Immunogen

Amino acids EPRKSDCAQDKS were used as the immunogen for this HOXB6 antibody.

Storage

Aliquot and store the HOXB6 antibody at -20oC.