


HOXB9 Antibody (R33217)

Catalog No.	Formulation	Size
R33217-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human, Mouse
Predicted Reactivity	Dog
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	3219
Applications	Western Blot : 0.3-1ug/ml ELISA (peptide) LOD : 1:4000
Limitations	This HOXB9 antibody is available for research use only.



Western blot testing of mouse spinal cord lysate with HOXB9 antibody at 0.3ug/ml.
Predicted molecular weight: ~28kDa.

Description

Additional name(s) for this target protein: Homeobox B9

Application Notes

Optimal dilution of the HOXB9 antibody should be determined by the researcher.

Immunogen

Amino acids CEGSEDKERPDQTN were used as the immunogen for this HOXB9 antibody.

Storage

Aliquot and store the HOXB9 antibody at -20oC.