


HOXC10 Antibody (R33944)

Catalog No.	Formulation	Size
R33944-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Predicted Reactivity	Pig, Cow
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	3226
Applications	Western Blot : 1-3ug/ml ELISA (peptide) LOD : 1:64000
Limitations	This HOXC10 antibody is available for research use only.



Western blot testing of HepG2 lysate with HOXC10 antibody at 1ug/ml. Expected molecular weight: 38-50 kDa.

Description

Additional name(s) for this target protein: Homeobox C10

Application Notes

Optimal dilution of the HOXC10 antibody should be determined by the researcher.

Immunogen

Amino acids EQRASLNPRAEHLE were used as the immunogen for this HOXC10 antibody.

Storage

Aliquot and store the HOXC10 antibody at -20oC.