


HOXC8 Antibody (R35973)


Catalog No.	Formulation	Size
R35973-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human, Mouse
Predicted Reactivity	Dog
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	3224
Applications	Western Blot : 0.1-0.3ug/ml ELISA (peptide) LOD : 1:2000
Limitations	This HOXC8 antibody is available for research use only.



Western blot testing of peripheral blood lymphocyte lysate with HOXC8 antibody at 0.1ug/ml. Expected molecular weight: 27-35 kDa.



Western blot testing of mouse spinal cord lysate with HOXC8 antibody at 0.1ug/ml. Expected molecular weight: 27-35 kDa.

Description

Additional name(s) for this target protein: Homeobox C8

Application Notes

Optimal dilution of the HOXC8 antibody should be determined by the researcher.

Immunogen

Amino acids SLYGAQQEAS were used as the immunogen for this HOXC8 antibody.

Storage

Aliquot and store the HOXC8 antibody at -20oC.