


## HOXC8 Antibody (R35973)


Catalog No.	Formulation	Size
R35973-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human, Mouse
Predicted Reactivity	Dog
Format	Antigen affinity purified
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	3224
Applications	Western Blot : 0.1-0.3ug/ml ELISA (peptide) LOD : 1:2000
Limitations	This HOXC8 antibody is available for research use only.



Western blot testing of peripheral blood lymphocyte lysate with HOXC8 antibody at 0.1ug/ml. Expected molecular weight: 27-35 kDa.



Western blot testing of mouse spinal cord lysate with HOXC8 antibody at 0.1ug/ml. Expected molecular weight: 27-35 kDa.

## Description

Additional name(s) for this target protein: Homeobox C8

## Application Notes

Optimal dilution of the HOXC8 antibody should be determined by the researcher.

## Immunogen

Amino acids SLYGAQQEAS were used as the immunogen for this HOXC8 antibody.

## Storage

Aliquot and store the HOXC8 antibody at -20oC.