


HOXD12 Antibody (R35166)

Catalog No.	Formulation	Size
R35166-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	3238
Applications	Western Blot : 2-3ug/ml ELISA (peptide) LOD : 1:1000
Limitations	This HOXD12 antibody is available for research use only.



Western blot testing of human heart lysate with HOXD12 antibody at 2ug/ml. Predicted molecular weight: ~29kDa, observed here at ~37kDa.

Description

Additional name(s) for this target protein: Homeobox D12

Application Notes

Optimal dilution of the HOXD12 antibody should be determined by the researcher.

Immunogen

Amino acids PTAKDGP^{EE}QAK were used as the immunogen for this HOXD12 antibody.

Storage

Aliquot and store the HOXD12 antibody at -20oC.