


HOXD12 Antibody (R35166)

Catalog No.	Formulation	Size
R35166-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Format	Antigen affinity purified
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	3238
Applications	Western Blot : 2-3ug/ml ELISA (peptide) LOD : 1:1000
Limitations	This HOXD12 antibody is available for research use only.



Western blot testing of human heart lysate with HOXD12 antibody at 2ug/ml. Predicted molecular weight: ~29kDa, observed here at ~37kDa.

Description

Additional name(s) for this target protein: Homeobox D12

Application Notes

Optimal dilution of the HOXD12 antibody should be determined by the researcher.

Immunogen

Amino acids PTAKDGPREEQAK were used as the immunogen for this HOXD12 antibody.

Storage

Aliquot and store the HOXD12 antibody at -20oC.