


HOXD13 Antibody (R34652)

Catalog No.	Formulation	Size
R34652-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Predicted Reactivity	Cow, Mouse, Rat
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	3239
Applications	Immunohistochemistry (FFPE) : 1-2ug/ml ELISA (peptide) LOD : 1:8000
Limitations	This HOXD13 antibody is available for research use only.



IHC testing of FFPE human prostate with HOXD13 antibody at 4ug/ml. HIER: steamed with pH6 citrate buffer, AP-staining.

Description

Additional name(s) for this target protein: Homeobox D13

Application Notes

Optimal dilution of the HOXD13 antibody should be determined by the researcher.

Immunogen

Amino acids KSSFPGDVALNQPD were used as the immunogen for this HOXD13 antibody.

Storage

Aliquot and store the HOXD13 antibody at -20oC.