


HPS3 Antibody (R35179)

Catalog No.	Formulation	Size
R35179-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Predicted Reactivity	Cow, Mouse, Pig, Rat
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	84343
Applications	Western Blot : 1-3ug/ml ELISA (peptide) LOD : 1:32000
Limitations	This HPS3 antibody is available for research use only.



Western blot testing of A431 lysate with HPS3 antibody at 1ug/ml. Predicted molecular weight: ~114kDa.

Description

Additional name(s) for this target protein: Hermansky-Pudlak syndrome 3, Cocoa

Application Notes

Optimal dilution of the HPS3 antibody should be determined by the researcher.

Immunogen

Amino acids PYLLYCSRKKPLT were used as the immunogen for this HPS3 antibody.

Storage

Aliquot and store the HPS3 antibody at -20oC.