


## HPS6 Antibody (R34739)

Catalog No.	Formulation	Size
R34739-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Predicted Reactivity	Cow, Mouse, Pig, Rat
Format	Antigen affinity purified
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	79803
Applications	Western Blot : 1-3ug/ml ELISA (peptide) LOD : 1:32000
Limitations	This HPS6 antibody is available for research use only.



Western blot testing of human lung lysate with HPS6 antibody at 1ug/ml. Predicted molecular weight: ~83kDa, observed here at ~95kDa.

## Description

Additional name(s) for this target protein: Hermansky-Pudlak syndrome 6

## Application Notes

Optimal dilution of the HPS6 antibody should be determined by the researcher.

## Immunogen

Amino acids LRTELRDYRGLEC were used as the immunogen for this HPS6 antibody.

## **Storage**

Aliquot and store the HPS6 antibody at -20oC.