


HRAS Antibody [clone BEA-8] (RQ4867)

Catalog No.	Formulation	Size
RQ4867	Antibody in PBS with 0.02% sodium azide, 50% glycerol and 0.4-0.5mg/ml BSA	100 ul

[Bulk quote request](#)

Availability	1-2 weeks
Species Reactivity	Human
Format	Purified
Host	Rabbit
Clonality	Rabbit Monoclonal
Isotype	Rabbit IgG
Clone Name	BEA-8
Purity	Affinity purified
UniProt	P01112
Applications	Western Blot : 1:500-1:2000
Limitations	This HRAS antibody is available for research use only.



Western blot testing of human MCF7 cell lysate with HRAS antibody. Predicted molecular weight ~21 kDa.

Description

HRAS is involved in the activation of Ras protein signal transduction.

Application Notes

Optimal dilution of the HRAS antibody should be determined by the researcher.

Immunogen

A synthetic peptide specific to human HRAS was used as the immunogen for the HRAS antibody.

Storage

Store the HRAS antibody at -20oC.