


Hras1 Antibody (F47429)

Catalog No.	Formulation	Size
F47429-0.4ML	In 1X PBS, pH 7.4, with 0.09% sodium azide	0.4 ml
F47429-0.08ML	In 1X PBS, pH 7.4, with 0.09% sodium azide	0.08 ml

Bulk quote request

Availability	1-3 business days
Species Reactivity	Mouse
Predicted Reactivity	Human, Rat, Chicken
Format	Antigen affinity purified
Clonality	Polyclonal (rabbit origin)
Isotype	Rabbit Ig
Purity	Antigen affinity
UniProt	Q61411
Applications	Western Blot : 1:1000
Limitations	This Hras1 antibody is available for research use only.



Hras1 antibody western blot analysis in mouse bladder tissue lysate. Predicted molecular weight ~ 21 kDa.

Description

Ras proteins bind GDP/GTP and possess intrinsic GTPase activity.

Application Notes

Titration of the Hras1 antibody may be required due to differences in protocols and secondary/substrate sensitivity.

Immunogen

A portion of amino acids 151-179 from the mouse protein was used as the immunogen for this Hras1 antibody.

Storage

Aliquot the Hras1 antibody and store frozen at -20oC or colder. Avoid repeated freeze-thaw cycles.