


HRASLS2 Antibody (R35247)

| Catalog No. | Formulation | Size |
|--------------|---|--------|
| R35247-100UG | 0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide | 100 ug |

[Bulk quote request](#)

| | |
|---------------------------|--|
| Availability | 1-3 business days |
| Species Reactivity | Human |
| Format | Antigen affinity purified |
| Host | Goat |
| Clonality | Polyclonal (goat origin) |
| Isotype | Goat Ig |
| Purity | Antigen affinity |
| Gene ID | 54979 |
| Applications | Immunohistochemistry (FFPE) : 5-10ug/ml ELISA (peptide) LOD : 1:16000 |
| Limitations | This HRASLS2 antibody is available for research use only. |



IHC testing of FFPE human liver with HRASLS2 antibody at 5ug/ml. HIER: microwaved with pH9 Tris/EDTA buffer, HRP-staining.

Description

Additional name(s) for this target protein: HRAS-like suppressor 2

Application Notes

Optimal dilution of the HRASLS2 antibody should be determined by the researcher.

Immunogen

Amino acids QELPYSLTSDNC were used as the immunogen for this HRASLS2 antibody.

Storage

Aliquot and store the HRASLS2 antibody at -20oC.