


# HRASLS2 Antibody (R35247)

Catalog No.	Formulation	Size
	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide $$	100 ug

### **Bulk quote request**

Availability	1-3 business days
Species Reactivity	Human
Format	Antigen affinity purified
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	54979
Applications	Immunohistochemistry (FFPE) : 5-10ug/ml ELISA (peptide) LOD : 1:16000
Limitations	This HRASLS2 antibody is available for research use only.



IHC testing of FFPE human liver with HRASLS2 antibody at 5ug/ml. HIER: microwaved with pH9 Tris/EDTA buffer, HRP-staining.

# **Description**

Additional name(s) for this target protein: HRAS-like suppressor 2

## **Application Notes**

Optimal dilution of the HRASLS2 antibody should be determined by the researcher.

### **Immunogen**

Amino acids QELPYSLTSDNC were used as the immunogen for this HRASLS2 antibody.

# **Storage** Aliquot and store the HRASLS2 antibody at -20oC.