


HRH1 Antibody (R33287)


Catalog No.	Formulation	Size
R33287-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Predicted Reactivity	Mouse, Rat, Dog, Cow
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	3269
Localization	Cytoplasmic
Applications	Immunohistochemistry (FFPE) : 4-6ug/ml ELISA (peptide) LOD : 1:32000
Limitations	This HRH1 antibody is available for research use only.



IHC testing of FFPE human colon with HRH1 antibody at 4ug/ml. HIER: steamed with pH6 citrate buffer, HRP-staining.



IHC testing of FFPE human small intestine with HRH1 antibody at 5ug/ml. HIER: steamed antigen retrieval with pH6 citrate buffer; AP-staining.

Description

Additional name(s) for this target protein: Histamine receptor H1

Application Notes

Optimal dilution of the HRH1 antibody should be determined by the researcher.

Immunogen

Amino acids CNENFKKTFKRILH were used as the immunogen for this HRH1 antibody.

Storage

Aliquot and store the HRH1 antibody at -20oC.