


## HRH3 Antibody (R36104)

Catalog No.	Formulation	Size
R36104-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Predicted Reactivity	Cow, Dog, Mouse, Rat
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
UniProt	Q9Y5N1
Gene ID	11255
Applications	Western Blot : 0.1-0.3ug/ml ELISA (peptide) LOD : 1:128000
Limitations	This HRH3 antibody is available for research use only.



Western blot testing of A549 lysate with HRH3 antibody at 0.1ug/ml. Predicted molecular weight: 36~49kDa (isoforms 1-7).

## Description

Additional name(s) for this target protein: Histamine H3 receptor

## Application Notes

Optimal dilution of the HRH3 antibody should be determined by the researcher.

## **Immunogen**

Amino acids SVTRAVSYRAQQGDT were used as the immunogen for this HRH3 antibody. This sequence is common to isoforms 1-7 (UniProt Q9Y5N1).

## **Storage**

Aliquot and store the HRH3 antibody at -20oC.