


HS3ST1 Antibody (R34842)

Catalog No.	Formulation	Size
R34842-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Predicted Reactivity	Cow
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	9957
Applications	Western Blot : 1-3ug/ml ELISA (peptide) LOD : 1:128000
Limitations	This HS3ST1 antibody is available for research use only.



Western blot testing of human cerebellum lysate with HS3ST1 antibody at 1ug/ml.
Predicted molecular weight: ~36kDa.

Description

Additional name(s) for this target protein: Heparan sulfate (glucosamine) 3-O-sulfotransferase 1

Application Notes

Optimal dilution of the HS3ST1 antibody should be determined by the researcher.

Immunogen

Amino acids NKLHEYFHEPNKK were used as the immunogen for this HS3ST1 antibody.

Storage

Aliquot and store the HS3ST1 antibody at -20oC.