


HSD11B1 Antibody (R33358)

| Catalog No. | Formulation | Size |
|--------------|-------------------------------------------------------------------------------|--------|
| R33358-100UG | 0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide | 100 ug |

[Bulk quote request](#)

| | |
|-----------------------------|-----------------------------------------------------------|
| Availability | 1-3 business days |
| Species Reactivity | Human |
| Predicted Reactivity | Cow |
| Format | Antigen affinity purified |
| Host | Goat |
| Clonality | Polyclonal (goat origin) |
| Isotype | Goat Ig |
| Purity | Antigen affinity |
| Gene ID | 3290 |
| Applications | Western Blot : 0.3-1ug/ml ELISA (peptide) LOD : 1:8000 |
| Limitations | This HSD11B1 antibody is available for research use only. |



Western blot testing of human liver lysate with HSD11B1 antibody at 0.3ug/ml. Predicted molecular weight: ~32kDa.

Description

Additional name(s) for this target protein: 11-beta-hydroxysteroid dehydrogenase 1

Researchers investigating glucocorticoid signaling, steroid metabolism, and cortisol activation pathways may also benefit from this [HSD11B1 Antibody / Glucocorticoid Metabolism Enzyme page](#) featuring knockout-validated and peptide-blocked endogenous tissue western blot data.

Application Notes

Optimal dilution of the HSD11B1 antibody should be determined by the researcher.

Immunogen

Amino acids CTSYNMDRFINK were used as the immunogen for this HSD11B1 antibody.

Storage

Aliquot and store the HSD11B1 antibody at -20oC.