


## HSD11B1 Antibody (R33358)

Catalog No.	Formulation	Size
R33358-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Predicted Reactivity	Cow
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	3290
Applications	Western Blot : 0.3-1ug/ml ELISA (peptide) LOD : 1:8000
Limitations	This HSD11B1 antibody is available for research use only.



Western blot testing of human liver lysate with HSD11B1 antibody at 0.3ug/ml. Predicted molecular weight: ~32kDa.

### Description

Additional name(s) for this target protein: 11-beta-hydroxysteroid dehydrogenase 1

### Application Notes

Optimal dilution of the HSD11B1 antibody should be determined by the researcher.

## **Immunogen**

Amino acids CTSYNMDRFINK were used as the immunogen for this HSD11B1 antibody.

## **Storage**

Aliquot and store the HSD11B1 antibody at -20oC.