


## HSD17B3 Antibody / 17-beta-HSD 3 (R35914)

Catalog No.	Formulation	Size
R35914-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Format	Antigen affinity purified
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	3293
Applications	Western Blot : 1-3ug/ml
Limitations	This HSD17B3 antibody is available for research use only.



Western blot testing of human testis lysate with HSD17B3 antibody at 2ug/ml. Predicted molecular weight: ~34kDa.

### Description

Additional name(s) for this target protein: 17-beta-hydroxysteroid dehydrogenase type 3

### Application Notes

Optimal dilution of the HSD17B3 antibody should be determined by the researcher.

### Immunogen

Amino acids SHFLNAPDEIQSLIH were used as the immunogen for this HSD17B3 antibody.

## Storage

Aliquot and store the HSD17B3 antibody at -20oC.