


HSH2D Antibody (R34287)

Catalog No.	Formulation	Size
R34287-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Predicted Reactivity	Mouse
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	84941
Applications	Western Blot : 0.3-1ug/ml ELISA (peptide) LOD : 1:4000
Limitations	This HSH2D antibody is available for research use only.



Western blot testing of K562 lysate with HSH2D antibody at 0.3ug/ml. Predicted molecular weight: ~39kDa.

Description

Additional name(s) for this target protein: Hematopoietic SH2 domain containing

Application Notes

Optimal dilution of the HSH2D antibody should be determined by the researcher.

Immunogen

Amino acids HSHVGYTLSYKAQ were used as the immunogen for this HSH2D antibody.

Storage

Aliquot and store the HSH2D antibody at -20oC.